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ENGLISH INSTRUCTIONS FOR USE

I. Introduction

KATANA Zirconia LT is a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia LT has 5 available thicknesses: (10mm, 14mm, 18mm, 22mm and 26mm) and has 3 shade variations: (LT10, LT12 and LT13). KATANA Zirconia LT is monochrome. KATANA Zirconia LT is recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.


II. Intended Use

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

III. Sintering Program


Sintering Program 1

Temperature	Programming Rate	Holding Time
Room Temp. – 1500°C (2732°F)	10°C/min. (18°F/min.)	—
1500°C (2732°F)	—	2 hrs
1500°C (2732°F) – Room Temp.	-10°C/min. (-18°F/min.)	—

Sintering Program 2

Temperature	Programming Rate	Holding Time
Room Temp. – 1515°C (2759°F)	35°C/min. (63°F/min.)	—
1515°C (2759°F)	—	0.5 hrs
1515°C (2759°F) – Room Temp. (*)	-45°C/min. (-81°F/min.)	—

IV. Composition

ZrO₃, Y₂O₃ etc.

V. Type and Class (ISO6872:2015)

Type: II/ Class: 5

VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 10.0×10⁻⁶/K

VII. Directions for Use

(1) Take the disc from the packaging and confirm that the disc does not have a crack or other damage.

(2) Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.

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

9. The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.

10. When preparing teeth, follow the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.

11. Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Anterior crown or bridge	0,4 mm or more
Veneer	0,4 mm or more
Posterior crown or bridge	0,5 mm or more
Inlay or onlay	0,5 mm or more

12. Use the following cross-sectional areas for connectors when fabricating bridges.

Location & indication	Connector cross section
Anterior	2- or 3-unit bridges more than 4-unit bridges
Posterior	2- or 3-unit bridges more than 4-unit bridges

13. After milling, remove the restorations from the disc with a diamond bur, etc.

(4) The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.

(5) Put the restorations into the refractory sagger tray and place them in the sintering furnace.

(6) Depending on the performance of the sintering furnace used, review the sintering schedule shown above (III. Sintering Program) before sintering the restorations.

(7) After sintering, adjust the restorations with a diamond bur as needed.

(8) Confirm that the restorations have no cracks.

(9)-1 FCZ:

a) Glaze baking: Create a high shine surface by polishing, especially on the contact areas, then apply the glaze on all surfaces in the usual manner.

b) Hand polishing: Create a high shine surface on the entire restoration by polishing without using the glaze.

(9)-2 Frameworks: Build up the dental porcelain (CERABIEN ZR or CZR PRESS LF) on the frameworks following the manufacturer's technical instructions. Check the coefficient of thermal expansion of the porcelain in the manufacturer's technical instructions to confirm compatibility.

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Caution:

1. This product should not be used when malocclusion, clenching or bruxism conditions are present.

2. When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.

3. When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.

4. Do not use for any purposes except for dental restoration. This product is for dental application only.

5. Do not touch the disc heated in the furnace with your bare hands.

6. Do not take the restorations out of the sintering furnace during high temperature, as the quenching causes the breaking.

However, if an elevating type furnace is used after sintering program 2, the restorations may be removed from the furnace at 800°C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.

7. There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2 is compatible.

8. Sintering program 2 is only recommended for frameworks (up to 3 units), FCZ crowns, FCZ bridges (up to 3 units), inlays, onlays and veneers.

[WARRANTY]

Kuray Noritake Dental Inc. will replace any product that is proven to be defective. Kuray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

"KATANA" and "CERABIEN" are trademarks of NORITAKE CO., LIMITED.

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